



A New Approach to Measure Contact Angle and Evaporation Rate With Flow Visualization in a Sessile Drop

Nengli Zhang
Ohio Aerospace Institute, Brook Park, Ohio

David F. Chao
Glenn Research Center, Cleveland, Ohio

Prepared for the
Interfaces for the Twenty-First Century: New Research Directions in Fluid
Mechanics and Materials Science
sponsored by the NASA OLMSA Microgravity Research Division, NSF
Monterey, California, August 16–18, 1999

National Aeronautics and
Space Administration

Glenn Research Center

Available from

NASA Center for Aerospace Information
7121 Standard Drive
Hanover, MD 21076
Price Code: A02

National Technical Information Service
5285 Port Royal Road
Springfield, VA 22100
Price Code: A02

A NEW APPROACH TO MEASURE CONTACT ANGLE AND EVAPORATION RATE WITH FLOW VISUALIZATION IN A SESSILE DROP

Nengli Zhang and David F. Chao
National Aeronautics and Space Administration
Glenn Research Center
Cleveland, Ohio 44135

The contact angle and the spreading process of sessile droplet are very crucial in many technological processes, such as painting and coating, material processing, film-cooling applications, lubrication, and boiling. Additionally, as it is well known that the surface free energy of polymers cannot be directly measured for their elastic and viscous restraints. The measurements of liquid contact angle on the polymer surfaces become extremely important to evaluate the surface free energy of polymers through indirect methods linked with the contact angle data.

Due to the occurrence of liquid evaporation is inevitable, the effects of evaporation on the contact angle and the spreading become very important for more complete understanding of these processes. Based on the laser shadowgraphic system used by the present author [1, 2], a very simple optical procedure has been developed to measure the contact angle, the spreading speed, the evaporation rate, and to visualize inside convection of a sessile drop simultaneously.

It is of interest to note that evaporation can induce Marangoni-Bénard convection in sessile drops [3]. However, the impacts of the inside convection on the wetting and spreading processes are not clear. The experimental methods used by previous investigators cannot simultaneously measure the spreading process and visualize the convection inside. The present work provides the solution for the first time.

The inside convection and the drop profile data can be synchronously recorded through the CCD camera I and II. CCD camera I would also precisely record the real-time diameter of the sessile drop which is essential for determination of both spreading speed and evaporation rate. As it is well known that a sessile drop takes the shape of a spherical cap provided its mass is sufficiently small, for example, less than 1 mg [4-6]. The sessile drop on the glass slide can be taken as a thin plano-convex lens with a focal length, f , which can be expressed as

$$f = \frac{R}{n - 1} \quad (1)$$

where R is the curvature radius of the spherical cap, n is the refractive index of the liquid. By a simple geometric relationship as shown in Fig. 1,

$$f = \frac{d(s + p)}{d + D} \quad (2)$$

where d is the base diameter of the sessile drop, s is the distance from the center of the drop base to its image on the front surface of the mirror, p is the distance between images of the drop center on the mirror and on the screen, D is the diameter of the drop image on the screen. Both d and D are time dependent due to the spreading and the evaporation, which can be accurately measured through the direct photo-shadowgraphic system. It is obvious that as the function of time, the curvature radius of the liquid spherical cap can be then determined by the data of d and D :

$$R(t) = (n - 1)(s + p) \frac{d(t)}{d(t) + D(t)} \quad (3)$$

The contact angle can be calculated as

$$\theta(t) = \arcsin \frac{d(t)}{2R(t)} \quad (4)$$

The evaporation rate can also be determined from the data of d and D . According to a simple geometric relationship, the instant volume of the sessile drop can be written as

$$V(t) = \pi h^2(t) \left[R(t) - \frac{h(t)}{3} \right] \quad (5)$$




Figure 1.—Schematic of experimental setup.

Therefore, the evaporation rate is

$$\dot{V} = 2\pi h \dot{h} \left(R - \frac{h}{3}\right) + \pi h^2 \left(\dot{R} - \frac{\dot{h}}{3}\right) \quad (6)$$

where

$$h = R - \sqrt{R^2 - d^2 / 4}, \quad (7)$$

$$\dot{h} = \dot{R} - \frac{R\dot{R} - d\dot{d} / 4}{\sqrt{R^2 - d^2 / 4}}, \quad (8)$$

and

$$\dot{R} = (n-1)(s+p) \left[\frac{\dot{d}}{d+D} - \frac{d(\dot{d}+\dot{D})}{(d+D)^2} \right] \quad (9)$$

As an example, Figs. 2 and 3 show the direct photography from top and the laser shadowgraphy at moment $t = 4$ sec., respectively, for an evaporating freon-113 sessile drop. Inside convection can be seen clearly in Fig.3.




Figure 2.—Direct photography from top. $d(t) = 5.72$ mm at $t = 4$ sec.




Figure 3.—Laser shadowgraphy. $d(t) = 85.7$ mm at $t = 4$ sec.

References

1. N. Zhang, and W.J. Yang, "Natural Convection in Evaporating Minute Drops," ASME Journal of Heat Transfer, vol. 104, pp. 656–662 (1982).
2. N. Zhang, and W.J. Yang, "Microstructure of Flow Inside Minute Drops Evaporating on a Surface," ASME Journal of Heat Transfer, vol. 105, pp. 908–910 (1983).
3. N. Zhang, and W.J. Yang, "Visualization of Evaporative Convection in Minute Drops by Laser Shadowgraphy," Rev. Sci. Instrum., vol. 54, no. 1, pp. 93–96 (1983).
4. R.G. Picknett, and R. Bexon, "The Evaporation of Sessile or Pendant Drops in Still Air," J. Colloid Interface Sci., vol. 61, no. 2, pp. 336–350 (1977).
5. S.M. Rowan, M.I. Newton, and G. McHale, "Evaporation of Microdroplets and the Wetting of Solid Surfaces," J. Phys. Chem., vol. 99, no. 35, pp. 13268–13271 (1995).
6. M.J. De Ruijter, J. De Coninck, and G. Oshanin, "Droplet Spreading: Partial Wetting Regime Revisited," Langmuir, vol. 15, no. 6, pp. 2209–2216 (1999).

REPORT DOCUMENTATION PAGE

Form Approved
OMB No. 0704-0188

Public reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to Washington Headquarters Services, Directorate for Information Operations and Reports, 1215 Jefferson Davis Highway, Suite 1204, Arlington, VA 22202-4302, and to the Office of Management and Budget, Paperwork Reduction Project (0704-0188), Washington, DC 20503.

1. AGENCY USE ONLY (<i>Leave blank</i>)		2. REPORT DATE December 1999	3. REPORT TYPE AND DATES COVERED Technical Memorandum	
4. TITLE AND SUBTITLE A New Approach to Measure Contact Angle and Evaporation Rate with Flow Visualization in a Sessile Drop			5. FUNDING NUMBERS WU-101-13-0A-00	
6. AUTHOR(S) Nengli Zhang and David F. Chao				
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES) National Aeronautics and Space Administration John H. Glenn Research Center at Lewis Field Cleveland, Ohio 44135-3191			8. PERFORMING ORGANIZATION REPORT NUMBER E-11966	
9. SPONSORING/MONITORING AGENCY NAME(S) AND ADDRESS(ES) National Aeronautics and Space Administration Washington, DC 20546-0001			10. SPONSORING/MONITORING AGENCY REPORT NUMBER NASA TM-1999-209636	
11. SUPPLEMENTARY NOTES Prepared for the Interfaces for the Twenty-First Century: New Research Directions in Fluid Mechanics and Materials Science sponsored by the NASA OLMSA Microgravity Research Division, NSF, Monterey, California, August 16-18, 1999. Nengli Zhang, Ohio Aerospace Institute, 22800 Cedar Point Road, Cleveland, Ohio 44142 (work funded under NASA Cooperative Agreement NCC3-620); David F. Chao, NASA Glenn Research Center. Responsible person, David F. Chao, organization code 6712, (216) 433-8320.				
12a. DISTRIBUTION/AVAILABILITY STATEMENT Unclassified - Unlimited Subject Categories: 25, 34, and 70 This publication is available from the NASA Center for AeroSpace Information, (301) 621-0390.			12b. DISTRIBUTION CODE Distribution: Nonstandard	
13. ABSTRACT (<i>Maximum 200 words</i>) The contact angle and the spreading process of sessile droplet are very crucial in many technological processes, such as painting and coating, material processing, film-cooling applications, lubrication, and boiling. Additionally, as it is well known that the surface free energy of polymers cannot be directly measured for their elastic and viscous restraints. The measurements of liquid contact angle on the polymer surfaces become extremely important to evaluate the surface free energy of polymers through indirect methods linked with the contact angle data. Due to the occurrence of liquid evaporation is inevitable, the effects of evaporation on the contact angle and the spreading become very important for more complete understanding of these processes. It is of interest to note that evaporation can induce Marangoni-Bénard convection in sessile drops. However, the impacts of the inside convection on the wetting and spreading processes are not clear. The experimental methods used by previous investigators cannot simultaneously measure the spreading process and visualize the convection inside. Based on the laser shadowgraphic system used by the present author, a very simple optical procedure has been developed to measure the contact angle, the spreading speed, the evaporation rate, and to visualize inside convection of a sessile drop simultaneously. Two CCD cameras were used to synchronously record the real-time diameter of the sessile drop, which is essential for determination of both spreading speed and evaporation rate, and the shadowgraphic image magnified by the sessile drop acting as a thin plano-convex lens. From the shadowgraph, the inside convection of the drop can be observed if any and the image outer diameter, which linked to the drop profile, can be measured. Simple equations have been derived to calculate the drop profile, including the instantaneous contact angle, height, and volume of the sessile drop, as well as the evaporation rate. The influence of the inside convection on the wetting and spreading processes can be figured out through comparison of the drop profiles with and without inside convection when the sessile drop is placed at different evaporation conditions.				
14. SUBJECT TERMS Spreading; Contact angle; Evaporation rate; Flow visualization			15. NUMBER OF PAGES 10	
			16. PRICE CODE A02	
17. SECURITY CLASSIFICATION OF REPORT Unclassified	18. SECURITY CLASSIFICATION OF THIS PAGE Unclassified	19. SECURITY CLASSIFICATION OF ABSTRACT Unclassified	20. LIMITATION OF ABSTRACT	